

GInger supports prediction of homologous recombination deficiency and patient response to PARPi treatment from shallow genomic profiles.

Pozzorini P, André G, Colletta T, et al. *Cell Rep Med.* 2023;4(12):101344. doi: 10.1016/j.xcrm.2023.101344

BACKGROUND



Homologous recombination deficiency (HRD) is a predictive biomarker for PARPi sensitivity in ovarian cancer.



In addition to BRCA1/2 testing, measuring genomic instability (GI) can expand identification of HRD-positive patients.




METHODS




RESULTS

Analytical performance in clinical samples (n=327)

High concordance




Low rejection rate



No significant difference in GI status OPA was observed between samples processed across **5 independent clinical laboratories, supporting a decentralized approach to HRD analysis.**

Clinical relevance: Retrospective analysis in subset of PAOLA-1 samples (n=195)


Patients classified by GInger™ as GI-positive had **significantly (over two times) longer progression-free survival (PFS)** when receiving olaparib + bevacizumab versus those receiving placebo + bevacizumab (hazard ratio 0.49; p=0.01).



GInger™ classification was non-inferior to the reference method for stratifying ovarian cancer patients with regards to PARPi treatment response.

CONCLUSIONS

- GInger™ is highly concordant with the centralized reference method for HRD assessment and shows lower rejection rate.
- GInger™ is an easy-to-implement deep learning method that allows accurate stratification of HRD samples based on lpWGS data.
- GInger™ can support the decentralized identification of ovarian cancer patients who may benefit from first-line maintenance treatment with PARPi.



ABBREVIATIONS

CI, confidence interval; FFPE, formalin-fixed paraffin-embedded; GI, genomic instability; GII, Genomic Integrity Index; HRD, homologous recombination deficiency; OPA, overall percentage agreement; PARPi, poly(ADP-ribose) polymerase 1 inhibitor; PFS, progression-free survival.

FURTHER READING

Buisson A, Saintigny P, Constantoulakis P, et al. Blinded-assessment of a solution to evaluate olaparib maintenance treatment efficacy in patients with ovarian cancer from the GINECO/ENGOT PAOLA-1 trial. *J Clin Oncol.* 2023;14:5588-5588. doi: 10.1200/JCO.2023.41.16_suppl.5588.

Pozzorini C, Andre G, Coletta C, et al. **GInger predicts homologous recombination deficiency and patient response to PARPi treatment from shallow genomic profiles.** *Cell Rep Med.* 2023 Dec 19;4(12):101344. doi: 10.1016/j.xcrm.2023.101344.

GInger™ data were generated using the SOPHiA DDM™ Dx HRD Solution, available as a CE-IVD product for In Vitro Diagnostic Use in Europe and Turkey only.

SOPHiA GENETICS products are for Research Use Only and not for use in diagnostic procedures unless specified otherwise.