



Enabling molecular laboratories to maximize the utility of Illumina's TST170 panel

SOPHiA DDM for TST170 is a fully integrated bioinformatic workflow (FASTQ to Report) for Illumina TruSight™ Tumor 170 panel. The SOPHiA DDM™ Platform combines analytical performance* with streamlined interpretation of complex genomic variants in the context of comprehensive genomic profiling. Its intuitive interface and capabilities have been developed based on the feedback from our 1000+ customers. It offers several features to facilitate variant assessment and maximize the utility of genomic detection.

Analytical performance*

Detection of:

- SNVs & Indels in 153 genes
- Partner-agnostic detection of 55 fusion genes with full coverage of NTRK 1,2,3 fusions and splice variants
- CNVs in 152 genes (versus 59 genes with standard analysis solutions)¹

Full access to analytical output:

- Ability to access source files and quality reports for analysis

The pipeline is compatible with Illumina NextSeq™ sequencers.

Analysis time from FASTQ files: 8 hours⁷

	Performance metrics ²	Observed	Lower CI (95%)
DNA Variants			
SNV & Indels	Sensitivity ³	100%	95.71%
	Precision ⁴	98.58%	97.81%
CNV/Amplifications ⁵	Sensitivity	100%	-
RNA Variants			
Fusions & Splice Variants ⁶	Sensitivity	96.38%	-
	True Positive	133/138	-
	False Negative	5/138	-

One Simple Intuitive Platform: Beyond Analytics

Accelerated assessment and reporting of genomic variants

The SOPHiA DDM™ Platform helps users to immediately focus on relevant genomic alterations. Several features facilitate their interpretation process:

- Hotspot screening
- Algorithm-supported variant pre-classification
- Fully customizable filters
- Comprehensive report

Confident decision-making

The OncoPortal Plus module for SOPHiA DDM™ matches tumor molecular profiles with clinical associations and available clinical trials, leveraging expertly curated evidence powered by JAX-CKB™. After interpretation, the flexible reporting tools enable users to prepare push-button, comprehensive reports that are customizable to their needs.

Product codes:

DL0734ILLRSM

Global support at every step

We offer local support available globally. Our Set Up Program provides assistance with assay setup for fast and worry-free transition to routine testing. All along the usage, our robust support including re-analysis and dedicated bioinformaticians help save time and resources, ensuring fast resolution of workflow disruptions.

Secure and unlimited data storage

The SOPHiA DDM™ Platform provides unlimited and unrestricted storage, while keeping data safe by applying the highest industrial standards of encryption in compliance with your local data security policies.

Access to SOPHiA GENETICS community

Through the SOPHiA DDM™ Platform genomics experts from >790 healthcare institutions interpret their findings and flag the pathogenicity level of variants. This highly valuable information enriches the variant knowledge base and is safely shared among the members of the community, supporting their decision-making process for research purposes.

Match genomic profiles with clinical trials opportunities

SOPHiA DDM™ Platform can support the identification of subjects who could benefit from investigational therapies. If the institution chooses to activate the feature, when a genomic profile matches a research opportunity, SOPHiA DDM™ Platform notifies relevant healthcare professionals. As a result, identified subjects are offered to participate in corresponding clinical trials, contingent on eligibility assessment and consent.

Somatic gene variant annotations and related content have been powered by, without limitation, The Jackson Laboratory Clinical Knowledgebase (JAX-CKB™).

***For Research Use Only. Not for Use in diagnostic procedures.**
OncoPortal Plus is for Clinical Decision Support Use Only. Not intended as a primary diagnostic tool.

All third party trademarks listed by SOPHiA GENETICS remain the property of their respective owners. Unless specifically identified as such, SOPHiA GENETICS' use of third party trademarks does not indicate any relationship, sponsorship, or endorsement between SOPHiA GENETICS and the owners of these trademarks. Any references by SOPHiA GENETICS to third party trademarks is to identify the corresponding third party goods and/or services and shall be considered nominative fair use under the trademark law.

1. Based on benchmarking studies comparing SOPHiA GENETICS with a standard analysis
2. Performance metrics calculated in samples processed on Illumina NextSeq™ 500/550
3. Based on analysis of 71 confirmed variants from 43 FFPE tumor samples from multiple tissues types
4. Based on analysis of 504 confirmed variants from AcroMetrix™ Oncology Hotspot Control
5. Based on analysis of 15 samples with confirmed gene amplifications
6. Based on analysis of 65 confirmed fusions from 105 clinical and reference samples
7. Analysis time may vary depending on the number of genes, samples analyzed and server load